



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

## 囊泡关联膜蛋白关联蛋白 A(VAPA)重组蛋白

适用生物 Homo sapiens (Human , 人)

产品应用 SDS-PAGE; WB; ELISA; IP. 规格 50ug

分子量 29.4kDa 价格 请咨询当地代理商

纯度 > 95%

来源 原核表达

## Recombinant Vesicle Associated Membrane Protein Associated Protein A (VAPA)

FOR IN VITRO AND RESEARCH USE ONLY, NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES!

Organism species	Homo sapiens (Human)
Product No.	YBF566Hu01
Source	Prokaryotic expression
Host	E. coli
Purity	> 95%
UOM	50ug
Predicted Molecular Mass	29.4kDa
Predicted isoelectric point	n/a
Applications	SDS-PAGE; WB; ELISA; IP.
Endotoxin Level	<1.0EU per 1μg (determined by the LAL method)

--

--

Subcellular Location	n/a
Residues	Met1~Leu249 (Accession # Q9P0L0) with N-terminal His-Tag
Formulation	Supplied as lyophilized form in PBS, pH7.4, containing 5% trehalose, 0.01% sarcosyl

## SEQUENCES



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

The target protein is fused with N-terminal His-Tag, its sequence is listed below.

MGHHHHHHSGSEF-MASASGAMAK HEQILVLDPP TDLKFKGPFT DVVTTNLKLR NPDRKVCFK V  
KTTAPRRYC VRPNSGIIDP GSTVTVSVML QPFDYDPNEK SKHKFMVQTI FAPPNTSDME AVWKE  
AKPDE LMDSKLRCVF EMPNENDKLN DMEPSKAVPL NASKQDGPMP KPHSVSLNDT ETRKLMEE  
CK RLQGEMMKLS EENRHLRDEG LRLRKVAHSD KPGSTSTASF RDNVTSPLPS LLVVIAAIFI GFF  
LGKFIL

## USAGE

Reconstitute in sterile PBS, pH7.2-pH7.4.

## STORAGE AND STABILITY

**Storage:** Avoid repeated freeze/thaw cycles. Store at 2-8°C for one month. Aliquot and store at -80°C for 12 months.

**Stability Test:** The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

## About the MARKER (complimentary)

**Effective Size Range:** 10kDa to 70kDa.

**Protein bands:** 10kDa, 14kDa, 18kDa, 22kDa, 26kDa, 33kDa, 44kDa and 70kDa.

**Double intensity bands:** The 26kDa, 18kDa, 10kDa bands are at double intensity to make location and size approximation of proteins of interest quick and easy.

**Ready-to-use:** No need to heat, dilute or add reducing agents before use.