

YBC395Hu01 100µg

Recombinant Carboxypeptidase A3, Mast Cell (CPA3) Organism Species: Homo sapiens (Human) Instruction manual

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

9th Edition (Revised in Jul, 2013)

[<u>PROPERTIES</u>]

kDa Residues: Ile110~Ser417 (Accession # P15088), with two 70 N-terminal Tags, His-tag and GST-tag. 44 Host: E. coli 33 26 Subcellular Location: Cytoplasmic vesicle, 22 secretory vesicle. 18 Purity: >95% 14 Endotoxin Level: <1.0EU per 1µg 10 (determined by the LAL method). Formulation: Supplied as lyophilized form in PBS, pH7.4 15% SDS-PAGE containing 5% trehalose, 0.01% sarcosyl. Predicted isoelectric point: 8.9 Predicted Molecular Mass: 68.4kDa Applications: SDS-PAGE; WB; ELISA; IP. (May be suitable for use in other assays to be determined by the end user.) [USAGE]



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

Reconstitute in sterile PBS, pH7.2-pH7.4.

[STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months. Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

[SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and GSTtag, its sequence is listed below.

MSPILGYWKI KGLVQPTRLL LEYLEEKYEE HLYERDEGDK WRNKKFELGL EFPNLPYYID GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV DFLSKLPEML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK KRIEAIPQID KYLKSSKYIA WPLOGWQATF GGGDHPPKSD GSTSGSGHHH HHHSAGLVPR GSTAIGMKET AAAKFERQHM DSPDLGTGGG SGIEGRGSMG YRGSEF-I PGRHSYAKYN NWEKIVAWTE KMMDKYPEMV SRIKIGSTVE DNPLYVLKIG EKNERRKAIF TDCGIHAREW VSPAFCQWFV YQATKTYGRN KIMTKLLDRM NFYILPVFNV DGYIWSWTKN RMWRKNRSKN QNSKCIGTDL NRNFNASWNS IPNTNDPCAD NYRGSAPESE KETKAVTNFI RSHLNEIKVY ITFHSYSOML LFPYGYTSKL PPNHEDLAKV AKIGTDVLST RYETRYIYGP IESTIYPISG SSLDWAYDLG IKHTFAFELR DKGKFGFLLP ESRIKPTCRE TMLAVKFIAK YILKHTS