TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com YB82576H

YB82576Hu01

# Fibronectin Type III Domain Containing Protein 5 (FNDC5)

Organism: Homo sapiens (Human)

Instruction manual


FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

5th Edition (Revised in January, 2013)

### [ DESCRIPTION ]

**Protein Names:** Fibronectin Type III Domain

#### **Human FNDC5**



Containing Protein 5 Synonyms: FNDC5

Species: Human

Size: 100µg

Source: Escherichia coli-derived

**Subcellular Location:** Cell membrane; Single-pass type I membrane protein. Peroxisome membrane;

Single-pass type I membrane protein.

## [ PROPERTIES ]

Residues: Leu2~Ile137 (Accession # Q8NAU1), with



TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

N-terminal His-Tag.

Grade & Purity: >95%, 19kDa as determined by

SDS-PAGE reducing conditions.

Formulation: Supplied as lyophilized form in PBS, pH

7.4, containing 5% sucrose, 0.01% sarcosyl.

Endotoxin Level: <1.0 EU per 1µg (determined by

the LAL method).

Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

Predicted Molecular Mass: 16.9kDa

Predicted isoelectric point: 6.2

TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

### [PREPARATION]

Reconstitute in sterile PBS, pH7.2-pH7.4.

### [STORAGE AND STABILITY]

Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

### [SEQUENCES]

The target protein is fused with N-terminal His-Tag, its sequence is listed below. MGHHHHHHSGSEF-**LRFIQEVNT TTRSCALWDL EEDTEYIVHV QAISIQGQSP** ASEPVLFKTP REAEKMASKN KDEVTMKEMG RNQQLRTGEV LIIVVVLFMW AGVIALFCRQ YDIIKDNEPN NNKEKTKSAS ETSTPEHQGG GLLRSKI