TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

YB92260Mu01 100µg Fusion (FUS) Organism: Mus musculus (Mouse) *Instruction manual* 

FOR IN VITRO USE AND RESEARCH USE ONLY NOT FOR USE IN CLINICAL DIAGNOSTIC PROCEDURES

7th Edition (Revised in May, 2013)


# [PROPERTIES]

Residues: Ser3~Arg267 (Accession # P56959), with two Nterminal Tags, His-tag and GST-tag. Host: *E. coli* Subcellular Location: Nucleus. Purity: >95% Endotoxin Level: <1.0EU per 1µg (determined by the LAL method). Formulation: Supplied as lyophilized form in PBS, pH7.4, containing 5% sucrose, 0.01% sarcosyl. Predicted isoelectric point: 7.7 Predicted Molecular Mass: 53.3kDa Applications: SDS-PAGE; WB; ELISA; IP.

(May be suitable for use in other assays to be determined by the end user.)

## [ <u>USAGE</u> ]

Reconstitute in sterile PBS, pH7.2-pH7.4.





TEL:4006-871-227 Web:www.ybio.net Email:shybio@126.com

## [STORAGE AND STABILITY]

#### Storage: Avoid repeated freeze/thaw cycles.

Store at 2-8°C for one month.

Aliquot and store at -80°C for 12 months.

Stability Test: The thermal stability is described by the loss rate of the target protein. The loss rate was determined by accelerated thermal degradation test, that is, incubate the protein at 37°C for 48h, and no obvious degradation and precipitation were observed. (Referring from China Biological Products Standard, which was calculated by the Arrhenius equation.) The loss of this protein is less than 5% within the expiration date under appropriate storage condition.

#### [SEQUENCES]

The target protein is fused with two N-terminal Tags, His-tag and GST-tag, its sequence is listed below.

MRNKKFELGL EFPNLPYYID GDVKLTQSMA IIRYIADKHN MLGGCPKERA EISMLEGAVL DIRYGVSRIA YSKDFETLKV DFLSKLPEML KMFEDRLCHK TYLNGDHVTH PDFMLYDALD VVLYMDPMCL DAFPKLVCFK KRIEAIPQID KYLKSSKYIA WPLQGWQATF GGGDHPPKSD **GSTSGSGHHHHHHSAGLVPRGSTAIGMKETAAAKFERQHMDSPDLGTLEV** LFQGPLGSEF- SNDYTQQA TQSYGAYPTQ PGQGYSQQSS QPYGQQSYSG YGQSADTSGY GQ S SY G S S Y G QTQ N T G Y G T Q SA P Q G Y G S T G GY G S S Q S S Q S SY G Q Q S S Y P G YGQ QPAPSST SGSYG GS SQS SSYGQ PQ SGG YGQ QSG YGG Q QQS YGQ QQ SS YNPPQGYGQQ NQYNSSSGGG GGGGGGGNYGQ DQSSMSGGGG GGGYGNQDQS GGGGGGYGGG QQDRGGRGRG GGGGYNRSSG GYEPRGRGGG RGGRGGMGGS DRGGFNKFGG PRDQGSR